

**QUICK DEFROSTING
EASY & CLEAN**

izidefrost+

DEFROSTING VACUUM

PROMOTE MAINTENANCE OF YOUR FREEZERS

-20°C / -80°C / -150°C

The essential tool for any laboratory with a fleet of freezers. Promote a more regular defrosting process, avoiding the need to turn off the freezer and transfer its contents. In doing so, extend the lifespan of your freezers and reduce your electricity consumption by **-30% à -50%**

● *Dimensions: L375mm x H53mm x W390mm - Weight: 9kg
Airflow: 4000l/m - Sound level: 64dba*

Professional device

Technical specifications specifically optimized for the aspiration of ice and snow at -80°C

Stainless steel scraper

Scrape and break the frost

25L tank

Vacuum the frost and ice

HEPA filter at the air outlet

Prevent any contamination in case of vacuuming a sample or sample residue



CONTACT@IZIBIO.FR

www.izibio.fr

